



The Universal Locking System

Only From Veterinary Implants Direct

Key Features

Patient Size				
Screw Insertion	All Screws Can Be Inserted in the Monoaxial or Polyaxial Orientation			
Screw Diameter	2.7 mm		3.5 mm	
Screw Type	Locking	ST Cortical	Locking	ST Cortical
Screw Driver	All Standard Hex Screwdrivers are Compatible with All Screws			
Drill Bit	2.0 mm	2.5 mm	3.2 mm	2.0 mm
Drill Guide	2.7 mm	3.5 mm	4.0 mm	2.7 mm
VID Locking Plates	Compatible with <u>All Screw Types and Sizes</u> (Locking, Non-Locking, 2.7mm, 3.5mm, 4.0mm)			
	 2.7 mm Limited Contact	 3.5 mm Locking Standard	 3.5 mm Locking Broad	
TPLO Locking Plates	Compatible with <u>All Screw Types and Sizes</u> (Locking, Non-Locking, 2.7mm, 3.5mm, 4.0mm)			
	 2.7 mm Universal Locking Plate	 3.5 mm 6-Hole Universal Locking Plate	 3.5 8-Hole Broad Universal Locking Plate	

Polyaxial Screw Insertion 15° in All Directions (30° Total)



No Need to Change Drill Bits Between Locking and Non-Locking Screws



Quick-Connect Drill Guide Design

Plates Accept All Screw Types and Sizes 2.7 mm - 4.0 mm



A Step Up in Efficiency

The autoclave-safe VID locking screw box is built from brushed aluminum and high temperature grade polymer resin. It houses 475 screws and includes a large instrument bay for plates, drill guides, and drill bits making your sterilization convenient and your time in the OR optimal.

www.VeterinaryImplantsDirect.com

